

Microlab[®] 600 Series GLP/GMP Compliance Package

Hamilton Company now offers an enhanced feature to the Microlab[®] 600 Advanced diluter and dispenser – a comprehensive compliance package. The diluter and dispenser can now be used in CFR 21 Part 11- and FDA GMP/ GLP-compliant laboratories, or any other facility that adheres to a security protocol, thanks to an accessory printer kit and new compliance software package that comes standard on all Advanced Microlab 600 devices and is available as a free software update for existing users. Highlights include the ability to administer user accounts and passwords, create log files that conform to CFR 21 Part 11 and manage log files on a PC using the LyncStore™ application.

The compliance assurance software update works in tandem with the new Microlab 600 Advanced printer kit. The kit is an accessory item that satisfies regulatory hardcopy requirements by allowing users to print log files directly from the device.

This package is designed for:

- ▶ FDA GXP regulated facilities
- ▶ Pharmaceutical manufacturing
- ▶ Clinical testing
- ▶ QA/QC laboratories
- ▶ Contract manufacturers



Compliance and Security Features



Enhanced Security Settings

- ▶ Create password protected user accounts, and define roles and permissions
- ▶ Track all user activities in log files compliant with CFR 21 Part 11
- ▶ Use electronic signatures and ensure user accountability



Advanced Log File Management

- ▶ Store log files on a 2GB card or permanently via LyncStore
- ▶ Archive, sort and filter log files with LyncStore PC-based management system
- ▶ Print log files to any printer connected to the PC



Convenient Printing

- ▶ Print log files directly from the controller
- ▶ Maintain internal records of pump usage
- ▶ Implement the requirements of close circle reporting

Ordering Information

Part Number	Description
61500-04	Printer Kit
ML625-DIL*	Dual Syringe Diluter with Advanced Controller
ML620-DIS*	Single Syringe Dispenser with Advanced Controller

* Note: Other models of the Microlab 600 are available



For more detailed information on the Microlab 600
GLP/GMP compliance package or to order, visit

www.hamiltoncompany.com/microlab600

© 2015 Hamilton Bonaduz AG. All rights reserved.
All trademarks are owned and/or registered by Hamilton Bonaduz AG.
Lit. No. 691087/02 – 12/2015

HAMILTON®

Web: www.hamiltoncompany.com

USA: 800-648-5950

Europe: +41-58-610-10-10

Hamilton Americas & Pacific Rim

4970 Energy Way
Reno, Nevada 89502 USA
Tel: +1-775-858-3000
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa

DJ 691 No. Km8+390m dreapta
307210 Giarmata/Romania
Tel: +40-356-635-050
Fax: +40-356-635-060
contact.lab.ro@hamilton.ch

To find a representative in your area, please visit hamiltoncompany.com.