



Ergonomic GL 45 cap

Retrace Code in cap

4 bottles per pack

Colourful Bottle Tags

Ergonomic bottle shape

Retrace Code on bottle

DURAN

Cap label

Nominal volume

Bottle label

Proven DURAN® glass

YOUTILITY – DESIGNED FOR YOU

THE NEW DURAN® LABORATORY BOTTLE SYSTEM

DURAN® YOUTILITY bottle GL 45, clear glass, complete

Cat. No.	Nominal Volume (ml)	Height (mm)	Outer diameter (mm)	Pack size
21 881 28 54	125	124	55	4
21 881 36 53	250	158	66	4
21 881 44 52	500	193	78	4
21 881 54 57	1000	253	93	4

DURAN® YOUTILITY bottle GL 45, amber glass, complete

Cat. No.	Nominal Volume (ml)	Height (mm)	Outer diameter (mm)	Pack size
21 886 28 59	125	124	55	4
21 886 36 58	250	158	66	4
21 886 44 57	500	193	78	4
21 886 54 53	1000	253	93	4

DURAN® YOUTILITY accessories and replacement parts

Cat. No.	Description	Material and colour	Pack size
29 229 28 02	DURAN® YOUTILITY Screw Cap GL 45	Food grade Polypropylene (PP), cyan	10
29 241 28 08	DURAN® YOUTILITY Pouring Ring GL 45	Food grade Polypropylene (PP), cyan	16
29 401 00 06	DURAN® YOUTILITY 200 self-adhesive label sets	White Polyester labels in a dispenser box	1 x 200 sets
29 243 29 04	DURAN® Bottle Tags GL 45 Mixed colours	Silicone. Eight colours: two bottle tags are supplied in each colour	16
29 243 28 18	NEU DURAN® Bottle Tags GL 45 Red	Silicone. Red colour only	20
29 243 28 26	NEU DURAN® Bottle Tags GL 45 Orange	Silicone. Orange colour only	20
29 243 28 34	NEU DURAN® Bottle Tags GL 45 Yellow	Silicone. Yellow colour only	20
29 243 28 42	NEU DURAN® Bottle Tags GL 45 Green	Silicone. Green colour only	20
29 243 28 59	NEU DURAN® Bottle Tags GL 45 Blue	Silicone. Blue colour only	20
29 243 28 67	NEU DURAN® Bottle Tags GL 45 Purple	Silicone. Purple colour only	20
29 243 28 75	NEU DURAN® Bottle Tags GL 45 Black	Silicone. Black colour only	20
29 243 28 83	NEU DURAN® Bottle Tags GL 45 White	Silicone. White colour only	20

For more information about the DURAN® YOUTILITY range get in touch at: sales@esslab.com